Microbial Cell Factories



Poster Presentation Open Access

Zea mays L. transglutaminase expression in Escherichia coli

Patricia Carvajal*¹, Enrique Villalobos², Alexandre Campos², J Ma Torné², Eduard Barberà¹ and Mireya Santos²

Address: ¹Department of Chemical Engineering, Biotechnology lab.-GQBB group, Institut Químic de Sarrià (IQS), Universitat Ramon Llull, Via Augusta 390, 08017, Barcelona, Spain and ²Institut de Biología Molecular de Barcelona-Consejo Superior de Investigaciones Científicas (IBMB-CSIC), Jordi Girona 18–26, 08034, Barcelona, Spain

* Corresponding author

from The 4th Recombinant Protein Production Meeting: a comparative view on host physiology Barcelona, Spain. 21–23 September 2006

Published: 10 October 2006

Microbial Cell Factories 2006, 5(Suppl 1):P8 doi:10.1186/1475-2859-5-S1-P8

© 2006 Carvajal et al; licensee BioMed Central Ltd.

Background

Transglutaminase (protein-glutamine:amine γ -glutamyltransferase, E.C. 2.3.2.13) catalyses acyl-transfer reactions between γ -carboxyamide groups of glutamine residues and the ϵ -amino group of lysines in proteins, leading to inter- or intramolecular cross-linking. Transglutaminases (TGs) have been found in mammals, plants, fish, nematodes and bacteria. Two maize cDNA clones (*TGZ15* and *TGZ21*) that expressed active transglutaminase localized in chloroplasts were isolated [1,2].

Results

A TGZ sequence was subcloned into pET (Novagen) vector (named as pET28a+*TGZ4*). The expression assays in *E. coli* BL21 DE3 cells showed that the main fraction of the protein (>80%) was found in the inclusion bodies (see Figure 1). The purification under denaturing conditions in FPLC system followed by a refolding step was a suitable procedure to obtain functional TGZ4p (see Figures 2,3). In addition, a specific antibody against TGZ4p was obtained in the laboratory for the immunolocalization of this protein in *E. coli* cells (see Figure 4).

Conclusion

TGZ4p was expressed in *E. coli* mainly as inclusion bodies. The purification under denaturing conditions and refolding *in vitro* was a suitable procedure to obtain functional TG.

Acknowledgements

Thanks to the Spanish National Project MCYT BFI2003-003318 and CEUB-IQS fellowship.

References

- Villalobos E, Santos M, Talavera D, Rodríguez-Falcón M, Torné J: Molecular cloning and characterization of a maize transglutaminase complementary DNA. Gene 2004, 336:93-104.
- Torné JM, Santos MA, Talavera D, Villalobos E: Maize nucleotide sequense coding for a protein with transglutaminase activity and use thereof. PCT/ES03/00247. Patent number WO03102128 2002.




Figure I
Western blot of TGZ4p. S = soluble, IB = inslusion bodies




Figure 2 TGZ4p FPLC purification under denaturing conditions.




Figure 3 TGase activity of refolded TGZ4p against Tris and phosphate buffer.




Figure 4 TGZ4p TEM Immunolocalization, S=soluble, IB=inclusion body.